

Technical Specifications



NOVEOS™ Dimensions

	Length	Height	Depth
NOVEOS™ Instrument	1562 mm (61.5 inches)	1295 mm (51 inches)	826 mm (32.5 inches)
Doors Open	2032 mm (80 inches)	1295 mm (51 inches)	1461 mm (57.5 inches)
Canopy Open		1867 mm (73.5 inches)	
Status Pole (not shown)		1625 mm (64 inches)	

NOTE: Allow an extra 914 mm (36 inches) to one side and to the front for servicing.



Weight

400 Kg (881 lbs.), with empty waste, wash and DI water containers.

Floor Space Requirement

Length	2489 mm (98 inches)
Depth	1651 mm (65 inches)

NOTE: Allow 305 mm (12 inches) at the back of the unit for power connection.

Desk/Bench Space Requirement

	Length	Depth	Height
Computer & Monitor Overall	914 mm (36 inches)	762 mm (30 inches)	409 mm (16.1 inches)
Computer	419 mm (16.5 inches)	419 mm (16.5 inches)	362 mm (14.25 inches)
Touchscreen Monitor LCD	521 mm (20.5 inches)	102 mm (4 inches)	409 mm (16.1 inches)
Keyboard	457 mm (18 inches)	152 mm (6 inches)	25 mm (1 inches)
Handheld Barcode Reader	66 mm (2.6 inches)	107 mm (4.2 inches)	168 mm (6.6 inches)

Power Requirements

Item	Requirements		
NOVEOS™ Instrument	100–240 VAC	50/60 Hz	1250W

Power Conditioner Power Conditioner

North America	North America	
Ametek Powervar Standard Conditioner	Input:	120VAC 60 Hz
ABC1200-11	Output:	120 VAC 60 Hz
	Input Plua:	5-15P

RoHS Compliant

Weight: 20 Kg (43 lbs.)

Output Receptacles: (6) 5-15R

Height: 283 mm (11.15 in.) Attention: Only NOVEOS™ System, PC and Touchscreen Monitor may be connected to the Power

Depth: 409 mm (16.1 in.) Conditioner

International International

Ametek Powervar Standard Conditioner Input: 230VAC 50/60 Hz
ABC1500-22INT Output: 230VAC 50/60 Hz

RoHS Compliant Input Plug: IEC 320
Output Receptacles: (9) IEC 320 C13

Weight: 20 Kg (43 lbs.)

Height: 283 mm (11.15 in.)

Width: 142 mm (5.6 in.)

Attention: Only NOVEOS™ System, PC and Touchscreen Monitor may be connected to the Power

Depth: 409 mm (16.1 in.) Conditioner.

60L for 1200 test run, when plumbed

Environmental Requirements

Operating Conditions Shipping and Storage Conditions Ambient Temperature Conditions Required without Fluids 18°C to 30°C (64.5°F to 86°F) **Ambient Temperature** Humidity 20 to 80% (non-condensing) -20°C to 50°C (-4°F to 122°F) relative humidity at 24°C (75°F) Humidity 20 to 95% (non-condensing) relative humidity at 24°C (75°F) Noise Output 20 dB **Heat Output** 4,265 BTU/hr. Waste Output 35L for 700 test run



Visit us online at

NOVEOSDx.com

HYCORbiomedical.com

These products are not currently available in the United States and in some regions outside of Europe.



