

Identify allergens with a system that gets it right

The right technology


- Quantitative total and specific IgE results with proven HYCOR enzyme immunoassay (EIA) system
- Comprehensive menu includes inhalant, food, occupational and environmental allergens as well as autoimmune testing
- Demonstrated correlation with fluorescence enzyme immunoassay (FEIA)¹
- Proven to support allergen immunotherapy in treatment of perennial allergic rhinitis and asthma²
- Used in major laboratories in U.S. and Europe

The right costs

- Get started easily and quickly and ensure smooth operation with outstanding HYCOR service
- Increase uptime and system utilization with proven HYTEC 288 Plus instrumentation
- Improve lab productivity and reduce labor costs with an automated, easy-to-use system

- ◀ HYCOR's proven activated cellulose solid-phase technology correlates with FEIA¹
- ▼ Windows-based software simplifies operation and enhances workflow management, data management and reporting of results

Ask your HYCOR representative for a cost analysis.



ANSWERS FOR ALLERGY CARE

HYCOR

An Agilent Technologies Division

Identify allergens with a system that gets it right

Get allergy testing right with HYCOR

Globally, HYCOR is partnering with leading labs like yours, providing high-quality products and setting high standards with quality customer service. Now, strengthened by the resources of Agilent Technologies, a premier measurement company, HYCOR is ready to partner with your lab to help you support your clients in the many facets of allergy care:

- Training and service: To help your lab get started and continue operating smoothly
- Outreach: A partnership with you to develop outreach programs, including physician education and patient communication materials for your physicians
- Website: Educational materials available through www.whatisallergy.com

hycorbiomedical.com

Agilent Technologies, Inc.

HYCOR Products Division
+1 800 382 2527 telephone
+1 714 933 3222 facsimile



HY•TEC™ 288 Plus: The system that gets it right

- Bench-top, solid-phase, automated enzyme immunoassay (EIA) instrument for quantitative total IgE and specific IgE
- Runs up to 50 patient samples, up to 288 results for each run
- Full allergen menu plus 27 autoimmune assays
- Manufactured in compliance with FDA's Quality System Regulation, current International ISO Standards, NCCLS guidelines; CE-marked

This information is subject to change without notice.

1. Kontis KJ, Valcour A, Patel A, et al. *Ann Clin Lab Sci.* 2006;36:79-87.
2. Yeoh KH, Wang DY, Gordon BR. *Otolaryngol Head Neck Surg.* 2004;131:673-678.

The HY•TEC 288 Plus instrument is designated for use by CLIA laboratories performing moderately complex assays
HYTEC is a trademark of Agilent Technologies, Inc.
© Agilent Technologies, Inc. 2009
Printed in USA January 30, 2009
5990-3498EN

ANSWERS FOR ALLERGY CARE

HYCOR

An Agilent Technologies Division